

rəversha by evoskin

LifePOWER EXOSOMES

The new era for Skin Regeneration

*From Cellular Longevity
to Visible Regeneration
Aesthetic Medicine
is entering a new era.*

***Skin Booster
for Face, Body & Scalp***



02

CAROMED ITALIA
MEDICAL RESEARCH

evoskin
YOU, AT YOUR BEST.

THE REGENERATIVE SKIN MANIFESTO

For decades, skin rejuvenation has relied on volume, inflammation, and mechanical stimulation. Fill, irritate, swell, repeat. Results were visible but often temporary, biologically superficial, and poorly aligned with true tissue health. We believe this approach is outdated. Real regeneration does not begin at the surface. It begins at the cell.

THE NEW PARADIGM: SEQUENCE MATTERS

Regeneration is not a single act. It is a biological sequence. We introduce a new regenerative logic

ACTIVATE LONGEVITY

Before asking the skin to regenerate, we must restore its responsiveness.

LongevHA™ is designed to:

- Reactivate SIRT1, the longevity gene
- Restore cellular energy through NAD⁺
- Reduce oxidative and inflammatory stress
- Rehydrate and stabilize the dermal matrix
- Prepare fibroblasts to respond intelligently to signals

This is not inflammation-based stimulation. This is cellular reprogramming. You cannot regenerate exhausted cells

DELIVER REGENERATION

Once the biology is restored, regeneration can be precise, efficient, and visible.

ReversHA™ is engineered to:

- Restore intercellular communication through certified milk exosomes
- Deliver PDRN, biomimetic peptides, and amino acids directly to fibroblasts
- Stimulate Type I and Type III collagen synthesis
- Modulate pigmentation pathways through glutathione and niacinamide
- Repair texture, tone, elasticity, and skin quality at the structural level

This is not volume. This is instruction-based regeneration. Cells don't heal without instructions.

WHY THIS IS DIFFERENT

This approach is: biological, not mechanical, signal-based, not trauma-based, compatible, not competitive, with devices and surgery, foundational, not temporary. It works with lasers, RF, microneedling, PRP, and surgery because it restores the biology those tools depend on.

The future of aesthetics is not about doing more. It is about doing things in the right order.

LongevHA™ turns the cell back on. ReversHA™ tells it what to do.

This is the new regenerative standard.

The power of exosomes



MILK EXOSOMES

- The most stable among any naturally occurring Exosomes
- Whitening & Antioxidant Effect
- Wound Healing
- Collagen Synthesis
- Anti-inflammatory Effect

2 MOLECULAR WEIGHT HA

- Keeps moisture & gives skin turgor

2 BIOMIMETIC OLIGOPEPTIDES

- Growth factor expression in skin cell lines
- Improves Hyperpigmentation

GABA PEPTIDE

- Stimulates expression
- Enhances production of HA

18 AMINO ACIDS

- All essential and almost all non essential ones
- Are necessary for maintaining healthy skin

PHARMACEUTICAL GRADE PRODUCTION

- ◆ **GMP Compliant Raw Materials**
- ◆ **Sterile Vials**
- ◆ **Endotoxin Tests for each batch**
- ◆ **Biocompatible**
- ◆ **Stable and Sterile Exosomes**



Protocol

Product Use Application Methods

- Cleanse the face thoroughly using a medical-grade disinfectant solution
- Apply up to 2.5 mL of the topical regenerative solution evenly over the treatment area.
- Perform microneedling using a pen device set at a 0.3 mm needle depth.
Alternatively, a transdermal ultrasonic device can be used to enhance penetration, applying the solution directly to the skin or via a cartridge.
- After microneedling, gently massage the area to ensure even distribution and absorption
- Perform four sessions spaced 2–3 weeks apart, based on clinical evaluation

Treatment Areas: Suitable for face, body, and scalp

Dosage: 2.5mL per zone

Post-Treatment: No social downtime; daily activities can resume immediately after the session

Standard protocol

- 3/4 sessions, one session every 3 weeks
- Maintenance: every 4 months at need

Product conservation

Store in a dry, cool and clean place at a temperature between 5°C and 25°C in the original packaging. Once opened, the product should be stored in the refrigerator and used within 14 days to ensure optimal efficacy and safety.

- **Certified Milk Exosomes (INCI)**
- **3 Year Shelf Life**
- **Skinbooster Effect**
- **Made in Italy**
- **360° Anti-age Formula**
- **Pharmaceutical Grade Production**
- **USA Origin Exosomes**
- **PDRN**



oversha



Clinical results that exceed expectations.

BEFORE



AFTER

BEFORE



AFTER

BEFORE



AFTER

Technical characteristics

MILK EXOSOMES		BOVINE
HIGH MOLECULAR WEIGHT & LOW MOLECULAR WEIGHT HYALURONIC ACID		10 MG/ML (2.4MDa-3MDa; 0.8-1MDa)
PDRN		5 MG/ML
PEPTIDES		GABA PEPTIDE OLIGOPEPTIDE-68 OLIGOPEPTIDE-24 & GLUTATHIONE
AMINO ACIDS	ALL ESSENTIAL & ARGININE GLYCINE PROLINE SERINE TYROSINE ALANINE ASPARTIC ACID GLUTAMIC ACID	
VITAMINS B GROUP		NIACINAMIDE PYRIDOXINE

REVERSHA Life Power Exosomes is a complete skinbooster for face, body and scalp produced according to pharmaceutical standards.

*Made in Italy
 5 vials for 5 mL
 1 pack for 10 treatments
 3 year shelf life*

TOPICAL USE



4300 Biscayne Blvd
 Miami FL 33137
 USA

T +1-802 386 7546
 E info@evoskin.com
 W evoskin.com

CAROMED ITALIA
 MEDICAL RESEARCH